


welfareplus+® Rodent Enrichment Solutions

SHELTER ENRICHMENT PRODUCT

Rabbit Shelter



Designed to provide a versatile hiding place that promotes well-being.

PRODUCT DESCRIPTION

The **Rabbit Shelter** provides a disposable, natural corrugated cardboard refuge, delivering a hygienic and stress-reducing sanctuary for primates and rabbits in research settings.

PRODUCT FEATURES AND BENEFITS

- **Stress Reduction & Social Housing:** Provides a secure and private space that helps reduce stress and fosters positive social interactions in group housing environments.
- **Ease of use:** Features a stackable design with an easy pop-up lock-in tab, allowing for quick and simple setup in any facility.
- **Guaranteed Hygiene:** Autoclavable material ensures thorough sterilization, maintaining a clean and safe environment for short-term studies.
- **Enhanced Activity & Visibility:** Multiple openings encourage natural exploratory behavior and allow for easy observation by researchers.
- **Ideal for Research Settings:** Suitable for short-term or toxicology studies and general research.
- **Suitable for ferrets and rabbits.**

PACKAGING

Product is autoclavable. Certificates of Analysis are available.

CODE	DESCRIPTION	QUANTITY
TE.5096C	Rabbit Shelter	50/bundle

Find out more

Contact your Inotiv representative
or visit us at [inotiv.com](https://www.inotiv.com)