

TD.130343

Control Rodent Diet

Formula

g/Kg

Casein	269.0
L-Cystine	4.0
Dextrose, monohydrate	444.235
Corn Oil	75.0
Lard	75.0
Cellulose	80.0
Mineral Mix, AIN-93G-MX (94046)	35.0
Vitamin Mix, AIN-93-VX (94047)	15.0
Choline Bitartrate	2.75
TBHQ, antioxidant	0.015

Footnote

Purified rodent diet containing ~27% casein, ~44% dextrose, and 7.5% each corn oil and lard. Although cited as a "fiber-free diet" (Desai et al., 2016, Cell 167, 1339–1353) it does contain 8% cellulose.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	23.8	23.6
Carbohydrate	42.7	42.4
Fat	15.3	34.0
Kcal/g	4.0	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotiv.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Dextrose
- Cellulose

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/tekla-orders
- teklaorders@inotiv.com
- 800.483.5523

Price Quotes | Order Status

- tekla@inotiv.com