

TD.03095

4% NaCl, Diet (Blue)

Formula	g/Kg
Wheat	350.0
Corn	263.79
Soybean Meal	197.0
Corn Gluten Meal (60%)	52.0
Alfalfa Meal (17%), dehydrated	30.0
Corn Oil	34.5
Sodium Chloride	40.0
Dicalcium Phosphate, FG (18.5% P, 21% Ca)	14.0
Calcium Carbonate, FG (38%)	12.0
Mineral Mix, TSD (80318)	1.5
Vitamin Mix, TSD (81125)	3.0
DL-Methionine, FG (99%)	1.0
L-Lysine HCl, FG (78%)	1.0
Ethoxyquin, antioxidant	0.01
Blue Food Color	0.2

Footnote

This is a modification of diet TD.92034 to add blue food coloring. This formula is part of a series of adjusted NaCl formulas of similar composition. Approximate nutrient content: 19% protein, 5% fat, 0.9% Ca, 0.65% P, 0.15% Mg, 0.8% K, 1.6% Na, and 2.5% Cl.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	18.5	24.2
Carbohydrate ²	45.9	60.1
Fat	5.4	15.8

Kcal/g 3.1

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com



Key Features

- Natural Ingredient Diet
- Adjusted NaCl Series
- Color Coded Blue
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com