

TD.170364

Doxycycline Diet (625, 2016, Red)

Formula

	g/Kg
2016, Teklad Global 16% Protein Rodent Diet	999.075
Doxycycline Hyclate	0.625
Red Food Color	0.3

Footnote

Diet adds 625 mg/kg doxycycline hyclate to Teklad's Global 2016 Rodent Diet and is color coded red. Diet was designed to deliver a daily dose of 2 - 3 mg of doxycycline based on consumption of 3 - 5 g/d by a mouse. Doxycycline hyclate contains approximately 87% doxycycline.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	16.2	21.1
Carbohydrate ²	52.8	68.6
Fat	3.5	10.4
Kcal/g	3.1	

¹ Values are calculated from ingredient analysis or manufacturer data

² Estimated digestible carbohydrate

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Standard Diet Base
- Doxycycline
- Teklad's Globa 2016
- Color Coded Red

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com