

TD.86489

Diet with Adjusted Sucrose/Corn Starch

Formula	g/Kg
Casein	207.0
DL-Methionine	3.0
Sucrose	320.99
Corn Starch	320.0
Corn Oil	50.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01

Footnote

A diet with approximately equal amounts (wt/wt) of sucrose and corn starch. About 53% of the carbohydrate kcal come from sucrose, and 47% from corn starch. The pellets are relatively short due to the 1:1 ratio of corn starch and sucrose.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	18.3	20.0
Carbohydrate	61.7	67.3
Fat	5.2	12.8
Kcal/g	3.7	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Adjusted Carbohydrate
- Sucrose & Corn Starch
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com