

TD.96355

Ketogenic Diet

Formula

g/Kg

Casein	173.3
DL-Methionine	2.6
Vegetable Shortening, hydrogenated (Crisco)	586.4
Corn Oil	86.2
Cellulose	87.97
Vitamin Mix, Teklad (40060)	13.0
Choline Bitartrate	2.5
TBHQ, antioxidant	0.13
Mineral Mix, Ca-P Deficient (79055)	20.0
Calcium Phosphate, dibasic	19.3
Calcium Carbonate	8.2
Magnesium Oxide	0.4

Footnote

A very high fat diet with almost no carbohydrate. The ratio of fat to protein + carbohydrate in this ketogenic diet is approximately 4.25.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	15.3	9.2
Carbohydrate	0.5	0.3
Fat	67.4	90.5
Kcal/g	6.7	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Ketogenic
- Crisco
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- Powder (paste)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com